



# Prism Coupler

---

Precise Measurement of  
Waveguide loss,  
Refractive Index  
and Film Thickness.

---

# SAIRON TECH

**The research of SaironTech (Established 1997) are specialized at the field industrial application and basic research for optical measurements.**

**SaironTech are developed(1998. 12) and commercialized the prism coupler (SPA- series) which is measuring the refractive index, thickness of the film and bulk material, and the waveguiding loss measurement.**



# Principle of Prism Coupler

Incidence Laser light

totally reflected at the base of the prism

Proper orientation of the direction of the incidence beam


Coupled through their evanescent fields in the gap.

Permits excitation of anyone of the film waveguide modes.

Measurement for both Refractive Index and Thickness of the film.

## Advantages of Prism Coupler

- **No advance knowledge required**
- **Unrelated Film/Substrate Combinations**
- **Dual-layer Film Measurement**
- **Bulk or Substrate Materials Measurement**
- **Anisotropy/Birefringence Measurement**
- **Thick- Film Measurement**



Make a narrow gab between the prism and the film  
using coupling (push) head

# Principle of Loss measurement

## Index Matching Oil Method

Waveguide light undergoes numerous total internal reflections inside the film.

Immerse the film into the liquid oil, with the index of refraction slightly higher than that of the film.


The light emerge out from the film at the intersection between the liquid surface and the film surface.

Detecting the outgoing light through the liquid oil.

Recording the intensity of the guided light as a function of propagation distance.


❖ **Highly precise measurement**

**down to 0.01dB/cm**



# Features of SPA-4000

- **Incidence laser light :: Rotation**
- **Prism and sample :: Fix**
- **Index matching oil method  
:: Loss measurement**
  - Sample size : Min. 1cm\*1cm to Max. 4-inch
  - Sample length for loss : more than 3cm
  - Loss measurable length : above 5cm



# Performance

Measurements	Specifications	
Refractive Index	Index measuring range	1.0 to 2.45
	Index accuracy	0.001
	Index resolution	$\pm 0.0005$
Thickness	Thickness measuring range	0.4 $\mu\text{m}$ ~ 20 $\mu\text{m}$
	Thickness accuracy	$\pm(0.5\%+50\text{\AA})$
	Thickness resolution	$\pm 0.3\%$
Bulk (index only)	Index accuracy	0.0005
	Index resolution	$\pm 0.0001$
Thick film (thickness only)	Thickness measuring range	2 $\mu\text{m}$ ~ 150 $\mu\text{m}$
Liquid (index only)	Index measuring range	1.0 to 2.4
	Index accuracy	$\pm 0.0005$
Loss Measurement	Measuring limitation	below 0.01dB/cm

# SPA-4000 Specification

SPA-series	include 632.8nm He-Ne Laser
	include GGG( $n=1.965$ ) prism & Holder ( index : $< 1.8$ )
	include One controller and PC interface(RS-232)
	include Analysis software (O/S : MS-windows_ Labview)
	include Si-photodiode Detector
Options	Laser Diode Module (405 ~1550nm); user choice
	TM Mode option for each wavelength
	Ge-Detector for INFRARED Laser
	Rutile( $n=2.865$ ) PRISM for high index refraction ( index : $1.8 \sim 2.45$ )
	THICK-film measurement system
	LIQUID measurement system
	Waveguide LOSS measurement system



# Application Fields

Film Types		Substrate Types
Silicon Nitride		Silicon
Silicon Dioxide	Silicon Oxynitride	GaAs
Low-k films		Quartz
Polymers	Polyimides	Glass
Zinc Sulfide	Titanium Dioxide	Sapphire
Sapphire	Epi Garnet	GGG
Photoresists	Holographic Gels	Lithium Niobate



# Application for SPA-4000

---

- **Optical components for optical communication systems**
  - **polymeric optical waveguide components**
    - **Optical Switches**
    - **Variable Optical Attenuators(VOA) for WDM(Wavelength Division Multiplexing)**
    - **Low optical propagation loss**
  - **Controllability of refractive index and birefringence**

**USERS : K-JIST( J.J. Kim, SPIE proceedings 7-11 July, Seattle, Washington, USA, 2002)**

<http://matlb.kjist.ac.kr/~optoelec/>

**ETRI ( M. H. Lee )** [http://www.etri.re.kr/e\\_etri/rnd/ard.html](http://www.etri.re.kr/e_etri/rnd/ard.html)

**Samsung**

<http://www.sait.samsung.co.kr/newsait/res/er.html>

**LG Electronics**

<http://www.lgelite.com/eng/aboutus/oe.html>

**Jeon-nam Univ.**

<http://physics.chonnam.ac.kr/bk21/>

- **Plastic Optical Fiber(POF)**
    - **Plastic Optical Fiber Amplifier(POFA) for Optical Communication**
    - **High Temperature Polymers for Waveguides**
- USERS : Hanyang Univ**  
[http://infochem.hanyang.ac.kr/index\\_intro.html](http://infochem.hanyang.ac.kr/index_intro.html)
-

# Application for SPA-4000

---

- **Properties of polymer**

- Investigation of chromic properties of polymer
- Information display and processing
- Storage Materials

**USERS** : KOREA AGENCY FOR TECHNICAL STANDARD (ATS)

[http://www.ats.go.kr/english/eng\\_home.asp](http://www.ats.go.kr/english/eng_home.asp)

**KRICT (E.K. Kim)**

[http://www.kRICT.re.kr/~cpl/index\\_e.html](http://www.kRICT.re.kr/~cpl/index_e.html)

**Hyosung**

<http://www.hyosung.co.kr/eng/index.jsp>

**LG Chemical**

[http://www.rnd.lgchem.co.kr/english/research/information/overview\\_set.htm](http://www.rnd.lgchem.co.kr/english/research/information/overview_set.htm)

- **Nano devices : MEMs, Micro-electronics**

**USERS** : ETRI

<http://std.etri.re.kr/eng/>

- **Temperature dependence**

**USERS** : K-JIST( J.J. Kim)



<http://matlb.kjist.ac.kr/~optoelec/>

---


# Application : Liquid index measurement

- \* Index range : 1.0 ~ 2.4
- \* Accuracy :  $\pm 0.0005$


- => Electromagnetic dependence
- => Temperature dependence
- => Birefringence



# Refractive index, Thickness Analysis GRAPH





# Waveguide LOSS Analysis GRAPH



# Layout for Prism Coupler

## Refractive Index, Thickness Measurement



## Loss Measurement

